DD Prefab

DE Gebrauchsanweisung, **EN** Instructions for use, **US** Instructions for use, **FR** Notice d'utilisation, **IT** Istruzioni per l'uso, **ES** Instrucciones de uso, **NL** Gebruiksaanwijzing, **EL** Οδηγίες χρήσης για το, **PT** Instruções de utilização, **TR** Kullanım talimatı, **RU** Инструкция по применению, **SV** Bruksanvisning, **ET** Kasutusjuhend, **DA** Brugsanvisning, **NO** Bruksanvisning, **FI** Käyttöohje, **IS** Notkunarleiðbeiningar, **PL** Instrukcja używania, **CS** Návod k použití skenovacího tělíska, **HU** Használati útmutató, **RO** Instrucțiuni de utilizare, **SK** Návod na použitie, **BG** Инструкция за употреба, **SR** Упутство за употребу, **FA** راهنمای استفاده از, **AR** دليل استخدام, **ZH** 使用说明书, **KO** 사용 설명서, **HR** Upute za uporabu, **GA** Treoracha úsáide, **LV** lietošanas pamācība, **LT** Naudojimo instrukcija, **SL** Navodila za uporabo, **MT** Struzzjonijiet għall-użu, **JA** 取扱説明書, **BS** Uputstvo za upotrebu

DE	DEUTSCH	Gebrauchsanweisung	4
EN	ENGLISH	Instructions for use	6
US	UNITED STATES (FDA)	Instructions for use	8
FR	FRANÇAIS	Notice d'utilisation	11
IT	ITALIANO	Istruzioni per l'uso	13
ES	ESPAÑOL	Instrucciones de uso	15
NL	NEDERLANDS	Gebruiksaanwijzing	17
EL	ΕΛΛΗΝΙΚΑ	Οδηγίες χρήσης για το	19
PT	PORTUGUÊS	Instruções de utilização	21
TR	TÜRKÇE	Kullanım talimatı	23
RU	РУССКИЙ	Инструкция по применению	25
SV	SVENSK	Bruksanvisning	27
ET	EESTI	Kasutusjuhend	29
DA	DANSK	Brugsanvisning	31
NO	NORSK	Bruksanvisning	33
FI	SUOMI	Käyttöohje	35
IS	ÍSLENSKA	Notkunarleiðbeiningar	37
PL	POLSKI	Instrukcja używania	39

CS	ČEŠTINA	Návod k použití skenovacího tělíska	41
HU	MAGYA	Használati útmutató	43
RO	ROMÂNĂ	Instrucțiuni de utilizare	45
SK	SLOVENSKY	Návod na použitie	47
BG	БЪЛГАРСКИ	Инструкция за употреба	49
SR	СРПСКИ	Упутство за употребу	51
FA	فارسی	راهنمای استفاده از	53
AR	عربي	دليل استخدام	55
ZH	中文	使用说明书	57
KO	한국어	사용 설명서	59
HR	HRVATSKI	Upute za uporabu	61
GA	GAEILGE	Treoracha úsáide	63
LV	LATVIEŠU VALODA	lietošanas pamācība	65
LT	LIETUVIŲ KALBA	Naudojimo instrukcija	67
SL	SLOVENŠČINA	Navodila za uporabo	69
MT	MALTI	Struzzjonijiet għall-użu	71
JA	日本語	取扱説明書	73
BS	BOSANSKI	Uputstvo za upotrebu	75

Inhalt:

- 1 x DD Prefab
- 1 x DD Lab Screw, blau
- 1 x DD Implant Screw, silberfarben

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- 1 x DD Prefab
- 1 x DD Lab Screw, blue
- 1 x DD Implant Screw, silver

Product description:

Dental Direkt Prefab blanks are prefabricated components for the machining of customized, one-piece abutments using CAD/CAM technology. They are available for various implant systems and consist of a prefabricated implant interface, screw channel, a retaining element for attachment in the milling machine and a groove for aligning the direction of rotation (axis of rotation). The cylindrical body located above the implant interface is transformed into the customized abutment design using computer-aided fabrication techniques. Dental Direkt Prefab blanks are made of Ti6Al4V (medical grade 5) for single use and are supplied non-sterile, each with a separately packaged Implant and Lab Screw.

Caution: The Dental Direkt Prefab blanks must only be used in combination with the DD Implant Library and the corresponding Prefab holder (incl. bar).

Caution: The Dental Direkt Prefab MDK blanks (DD Prefab MDK) must only be used in combination with the original Medentika® Implant Library and the corresponding original Medentika® Prefab holder.

Intended purpose:

DD Solid Connect® implant components are used to support prosthetic restorations in combination with endosseous dental implants in the upper and/or lower jaw.

Indication:

For fabrication of individual abutments on implants in the upper and lower jaw.

Contraindication:

- Intolerance to the constituents
- Restorations with an angulation correction of more than 30° to the implant axis

Please note: Observe the contraindications specified by the original implant manufacturer used in each case.

Intended users:

DD Solid Connect® implant components should only be used by dentists and dental technicians in compliance with the procedures in the instructions for use.

Intended patient group:

Fixed restorations using Dental Direkt Solid Connect® implant components are suitable for the permanent dentition in adult patients of any gender and nationality.

Instructions for use in the laboratory:

Make sure that the holder approved by Dental Direkt (DD Prefab holder and bar) or the original Medentika® Prefab holder for the Dental Direkt Prefab MDK blanks (DD Prefab MDK) is used to clamp the blank in the CAD/CAM system. This holder ensures that the implant sits in the machine vibration-free and is precisely orientated in the insertion direction. The prefabricated implant interface must not be modified during further processing. Damaged Prefabs must not be used. Always use suitable cooling and an extraction system when machining the blanks. Metal dusts are harmful to health, so always wear respiratory protection and safety goggles.

The Prefab blank must not show any signs of damage. An undamaged and correct implant interface is essential for correct placement of the abutment.

Caution: When designing in CAD, ensure that the minimum wall thickness (0.4 mm) is observed.

Caution: If not machined properly, titanium can overheat and ignite.

Caution: Only use the Lab Screw in the model (max. torque: 10 Ncm).

Instructions for use in the dental practice:

The individual abutment is inserted with the silver Implant Screw into the cleaned endosseous dental implant located in the patient's mouth.

Caution: The Lab Screw must not be used for placement of the implant components in the patient's mouth.

Caution: The Implant Screw must be tightened to the torque specified for the system (see table).

Please note: It is recommended to tighten the Implant Screw 10 minutes after placement.

Please note: When DD Solid Connect® implant components are used for restorations with reduced-diameter implants, limited indications apply. These can be found in the relevant instructions for use from the implant manufacturer.

Please note: DD Solid Connect® implant components from a particular series can only be combined with the corresponding compatible implant system. This also applies to the corresponding platform diameters.

Material:

DD Prefab, DD Implant Screw and DD Lab Screw are all made of Ti6Al4V (medical grade 5) ASTM 136 ELI.

Cleaning, disinfection and sterilization:

DD Prefabs and accompanying screws are cleaned at Dental Direkt, but not sterilized, before being packaged. After machining, and before being used on the patient, the implant components must be cleaned, disinfected and sterilized following the procedure described by Dental Direkt.

Cleaning/disinfection:

Automatic cleaning in a suitable washer-disinfector with an appropriate cleaning program for medical instruments (e.g. "Vario TD Program", Miele).

OR

The following manual cleaning procedure:

- Rinse with cold tap water for 1 minute
- Rinse the screw channel using a disposable syringe (min. capacity 10 ml)
- Clean in the ultrasonic cleaner at 60°C (140°F) +/- 2°C, e.g. using Dr. Weigert neodisher® LM 2 for 10 minutes at a concentration of 5%
- Neutralize with cold distilled water for 2 minutes
- Rinse in the ultrasonic cleaner using distilled water at 60°C (140°F) for 10 minutes

Sterilization:

Double-seal products in sterilization film conforming to EN ISO 11607-1. Steam sterilize in a suitable sterilizer with 3x fractionated prevacuum: 134°C (273°F) for 5 minutes, drying time: 20 minutes.

Then store in sterile packaging.

Caution: DD Prefab, DD Implant Screw and DD Lab Screw are suitable for single use only. Any reuse may lead to hazards due to infection.

Possible side effects and interactions:

Allergic reactions or skin irritations are possible in very rare cases. Different alloys in the oral cavity can lead to galvanic reactions when in direct contact.

Disposal:

In compliance with local regulations.

Warranty:

Our products undergo strict quality control and conform to the latest technical standards. Please refer to our current warranty terms for further details.

Please note: Any serious incident that occurs in connection with the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located

Please note: The safety and clinical performance summary report can be requested at info@dentaldirekt.de.

We are continuously developing and enhancing our devices, and therefore reserve the right to make changes. The latest version of the instructions for use can also be found on our website at: www.dentaldirekt.de/en

This version replaces all previous versions.

Explanation of symbols:

Manufacturer

Date of manufacture

Batch code

Catalog number

Keep dry

Non-sterile

Do not use if package is damaged

Torque

Consult electronic instructions for use
www.dentaldirekt.de/IFU

Caution: Under U.S. federal law, the device may only be sold by or on behalf of a dentist.

Do not re-use

Medical device

Unique device identification

Dental Direkt Series	Compatible with: manufacturer	System	Size (mm)	Torque (Ncm)
DD AAS-Serie	DENTSPLY Implants®	Astra Tech®	3.5 / 4.0	25
			4.5 / 5.0	25
DD BSK-Serie	Bredent Medical	Sky®	3.5 – 5.5	25
DD CCA-Serie	Altatec	Camlog®	3.3	20
			3.8	20
			4.3	20
			5.0	20
DD CCO-Serie	Altatec	Conelog®	3.3	20
			3.8	20
			4.3	20
			5.0	20
DD DEV-Serie	DENTSPLY Implants®	Astra EV	3.0	25
			3.6	25
			4.2	25
			4.8	25
			5.4	25
DD DFX-Serie	DENTSPLY Implants®	Frialit®/Xive®	3.4	25
			3.8	25
			4.5	25
			5.5	25
DD MAR-Serie	MegaGen	AnyRidge®	3.5 – 8.0	30
DD MMI-Serie	Medentis Medical®	ICX	3.75 – 4.8	30
DD MSE-Serie	MIS Implants®	Seven®	3.3 NP	30
			3.75 / 4.2 SP	30
			5.0 WP	30
DD NBR-Serie	Nobel Biocare™	Brånemark™	3.5 NP	35
			4.1 RP	35
DD NNA-Serie	Nobel Biocare™	Nobel Active®	3.5 NP	25
			4.3 / 5.0 RP	35
DD NRE-Serie	Nobel Biocare™	Replace Select®	3.5 NP	35
			4.3 RP	35
			5.0 WP	35
DD OTS-Serie	Osstem Implant	TS	3.5 M	20
			4.0 – 6.0 R	20
DD SBO-Serie	Straumann®	Bone Level®	3.3 NC	35
			4.1 / 4.8 RC	35
DD SSY-Serie	Straumann®	Tissue Level	4.8 RN	35
			6.5 WN	35
DD TSP-Serie	Thommen Medical	SPI®	4.0	25
			4.5	25
			5.0	25
			6.0	25
DD ZBC-Serie	Zimmer Biomet 3i	Certain®	3.4	30
			4.1 / 5.0	30
DD ZTA-Serie	Zimmer Dental	Tapered Screw-Vent®	3.5	30
			4.5	30
			5.7	30

Content:

- 1 x DD Prefab
- 1 x DD Lab Screw blue
- 1 x DD Implant Screw, silver

Indications for Use:

DD Solid Connect® CAD/CAM Abutments are used to support prosthetic restorations in combination with endosseous dental implants in the upper and/or lower jaw.

All digitally designed custom abutments for use with DD Solid Connect® CAD/CAM Abutments are to be sent to a Dental Direkt validated milling center* for manufacture.

DD Prefab abutments, for the Zimmer Biomet 3i Certain 3.4mm implant bodies, are indicated for maxillary lateral and mandibular central/lateral incisors only.

Compatible Implant Systems:

Dental Direkt Series	Compatible with: manufacturer	System	Size (mm)	Torque (Ncm)
DD AAS-Serie	DENTSPLY Implants®	Astra Tech®	3.5 / 4.0	25
			4.5 / 5.0	25
DD CCA-Serie	Altatec	Camlog®	3.3	20
			3.8	20
			4.3	20
			5.0	20
DD DEV-Serie	DENTSPLY Implants®	Astra EV	3.0	25
			3.6	25
			4.2	25
			4.8	25
			5.4	25
DD DFX-Serie	DENTSPLY Implants®	Frialit®/Xive®	3.4	25
			3.8	25
			4.5	25
			5.5	25
DD NNA-Serie	Nobel Biocare™	Nobel Active®	3.5 NP	25
			4.3 / 5.0 RP	35
DD NRE-Serie	Nobel Biocare™	Replace Select®	3.5 NP	35
			4.3 RP	35
			5.0 WP	35
DD SBO-Serie	Straumann®	Bone Level®	3.3 NC	35
			4.1 / 4.8 RC	35
DD SSY-Serie	Straumann®	Tissue Level	4.8 RN	35
			6.5 WN	35
DD ZBC-Serie	Zimmer Biomet 3i	Certain®	3.4	30
			4.1 / 5.0	30
DD ZTA-Serie	Zimmer Dental	Tapered Screw-Vent®	3.5	30
			4.5	30
			5.7	30

Device description:

The Dental Direkt Prefab blanks are prefabricated components for the metal-cutting production of individualized and one-piece abutments using CAD/CAM technology. The Dental Direkt Prefab blanks are available for various implant systems and consist of a prefabricated implant connection (interface), screw channel, a retaining element for fastening in the milling machine and a groove for aligning the direction of rotation (axis of rotation). The cylindrical body located above the implant connection is transformed into the individualized abutment design using computer-aided manufacturing techniques. The Dental Direkt Prefab blanks are made of Ti6Al4V (medical grade 5) for single use and are supplied non-sterile, each with a separately packed implant and laboratory screw.

Caution: The Dental Direkt prefab blanks must only be used in conjunction with the DD Implant Library and the corresponding prefab holder (incl. bar).

Caution: The Dental Direkt Prefab MDK blanks (DD Prefab MDK) must only be used in combination with the original Medentika® Implant Library and the corresponding original Medentika® Prefab holder.

Warning: Small diameter implants and angled abutment are recommended for the anterior region.

Indication:

DD Prefab: For fabricating customized abutments on implants in the upper and lower jaw.

DD Implant Screw: The implant screw (abutment screw) is used to fix the abutment to the implant (two-part implant system).

Material:

The DD Prefab, DD Lab Screw and DD Implant Screw are made of Ti6Al4V (medical grade 5) ASTM 136 ELI.

Contraindication:

- Intolerance to the components contained
- Restorations with an angulation correction of more than 20° to the implant axis

Intended users:

The DD Solid Connect® Implant Components may only be processed by trained, dental technical / dental medical personnel in compliance with the specifications in the Instructions for Use.

Intended patient group:

Fixed dentures made of Dental Direkt Solid Connect® CAD/CAM Abutments are suitable for the permanent dentition of adult patients of any gender and nationality.

Possible side effects and interactions:

Allergic reactions or skin irritations are possible in very rare cases. Different alloys in the oral cavity can lead to galvanic reactions when in direct contact.

Use in MR (Magnetic Resonance) environment:

MR Conditional

Warning: The RF Safety of the device has not been tested. The patient may only be imaged by landmarking at least 30 cm from the implant, or ensuring the implant is located outside the RF coil.

A patient with this device can be scanned safely in an MR system under the following conditions:

Device Name	DD Prefab, DD Implant screw
Static Magnetic Field Strength (B0)	≤ 3.0T
Maximum Spatial Field Gradient	30 T/m (3,000 gauss/cm)
RF Excitation	Circularly Polarized (CP)
RF Transmit Coil Type	For body transmit coil, landmarking at least 30 cm from the implant, or ensuring the implant is located outside of the coil. Extremity T/R coils permitted. Excludes Head T/R coil.
Operating Mode	Normal Operating Mode in the allowed imaging zone
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	Not evaluated for head landmark
Scan Duration	No specific constraints due to implant heating

Cleaning, disinfecting and sterilizing:

The DD Prefabs and the associated implant screws are cleaned at Dental Direkt, but packed non-sterile. After processing and before use on the patient, the implant components must be cleaned, disinfected and sterilized following the procedure described by Dental Direkt.

Cleaning / disinfecting:

Automatic cleaning in a washer-disinfector suitable for this application with an appropriate cleaning program for medical instruments (e.g. "Vario TD Program", Miele).

The following steps and parameters are recommended:

- Pre-washing step with cold tap water for 1 minute
- Washing step at 55°C +/- 1°C (131°F +/- 2°F) for 10 minutes, e.g. using Dr. Weigert neodisher® MediClean Dental
- Rinse twice (2 separate steps):
 1. First rinsing for 1 minute with cold deionized water <40°C (<104°F)
 2. Second rinsing for 2 minutes with cold deionized water <40°C (<104°F)
- Thermal disinfection step with 93°C (200°F) for 5 minutes
- Drying step at 110°C +/- 2°C (230°F +/- 4°F) for 20 minutes

OR

the following manual cleaning procedure:

- Rinse with cold tap water for 1 minute
- Rinse the screw channel using a disposable syringe (min. capacity 10 ml)
- Clean at 60°C +/- 2°C, e.g. using Dr. Weigert neodisher® LM 2 for 10 minutes at a concentration of 5% in the ultrasonic cleaner
- Neutralize with cold distilled water for 2 minutes
- Rinse in the ultrasonic cleaner using distilled water at 60°C for 10 minutes

Sterilizing:

Double-seal the product in FDA cleared sterilization film conforming to ISO 11607-1 and sterilize according to validated parameters listed below.

Steam sterilize in a FDA cleared sterilizer with 3x fractionated prevacuum: 132°C (270°F) for 4 minutes, drying time: 20 minutes.

Caution: After sterilization, the product must be used immediately and must not be stored. DD Prefab, DD Implant Screw and DD Lab Screw are suitable for single use only. Reuse may result in hazards due to infections.

Workflow:

1. Conventional impression or intra-oral scanning

Impression of the teeth setup could be done by conventional impression technique or by intra-oral scanning.

DD Scanbody 2.0 should be used for both.

The DD Scanbody 2.0 is used to record the position of the implant in the mouth, or on the model during the scanning process. Both the position of the implant and the rotation are depicted due to the asymmetric form of the DD Scanbody 2.0. It is important to select the correct DD Implant Library for referencing the position.

The DD Scanbody 2.0 is supplied for use on the model with a laboratory screw and with a silver-coloured implant screw for use in the patient's mouth. The surface must be clean and free of grease before scanning.

1.1 Conventional impression (Semi-digital workflow)

Take a conventional impression of the patient's tooth set-up and create a working model with the analog inserted to represent the implant.

Place the DD Scanbody 2.0 into the analog (i.e. DD Fusion analog) to determine the position and orientation of the implant.

Scan the working model by use of a desktop scanner (i.e. MEDIT T710).

1.2 Intra-oral scanning (Digital workflow)

To detect the exact implant position, use the DD Scanbody 2.0.

Selection must be carried out according to the compatibility of the implant system and the platform diameter.

Scan the patient's teeth setup by use of a dental intra-oral scanner (i.e. MEDIT i700).

2. Designing

The customized abutment must be designed using CAD Software Dental Direkt designer by exocad® (DentalCAD; K193352) with CAD Module DD Abutment Creator (Implant Module by exocad®, K193352) with the relevant library file installed.

CAD Software and CAD Module can be purchased from our website: www.dental-direkt.de/en

The library file can be downloaded from: www.dentaldirekt.de/en/downloads

Operation manual and tutorial videos for exocad® DentalCAD and Implant Module can be found on: www.wiki.exocad.com/wiki/index.php/Main_Page

Import the digitalized patient data to the design software (CAD Software Dental

Direkt designer by exocad® with CAD Module DD Abutment Creator). It is recommended to use always the latest version of the software.

Use appropriate library file and select relevant implant platform from the library.

Design the abutment in the design software according to the parameters.

The following parameters must be observed for the customized abutment:

Parameter	Specification
Abutment Post Height	4 mm – 12.5 mm
Margin height	0.5 mm – 6 mm
Diameter	12 mm maximum
Wall thickness	0.5 mm minimum
Angle from axis of implant	0° - 20°

3. Milling

All digitally designed custom abutments for use with DD Solid Connect® CAD/CAM Abutments are to be sent to a Dental Direkt validated milling center* for manufacture.

Make sure that the holder approved by Dental Direkt (DD Prefab holder and bar) or the original Medentika® Prefab holder for the Dental Direkt Prefab MDK blanks (DD Prefab MDK) is used to clamp the blank in the CAD/CAM system. This holder ensures the blank sits in the machine without vibrating, and the insert direction of the implant is oriented precisely. The prefabricated implant interface may not be amended during subsequent processing. Nor may damaged prefabs be used. Processing blanks in the machine is only approved with use of appropriate cooling. After milling, the abutment will be sent back to the dental lab for finalization.

Caution: Titanium can overheat and ignite if not processed properly.

4. Finalizing

Place the DD Prefab with the DD Lab Screw (blue) on the working model. Complete the crown/bridge restoration following routine laboratory procedures.

Caution: For use in the model, only use the laboratory screw (max. torque: 10 Ncm).

It is imperative to ensure adequate suction when processing.

Metal dust is harmful to the health, so always use breathing protection and protective glasses.

The prefab blank must not show any damage. Only an intact and correct implant interface ensures correct use of the abutment.

5. Use and handling for the dentist

The dentist receives the final dental restoration / working model with the DD Lab Screw (blue) from the dental lab.

Clean, disinfect and sterilize the dental restoration according to the instructions above.

Insert the dental restoration into the patient's mouth with the silver-coloured implant screw (DD Implant Screw) using a dental torque wrench with a suitable screw driver.

The implant screw should be tightened with the system - specific torque.


The screw channel must always be sealed after the abutment is attached to the implant.

Caution: The lab screw must not be used for inserting the implant components in the patient.

Note: Restricted indications apply for the restoration of DD Solid Connect® CAD/CAM Abutments on reduced-diameter implants, which can be found in the corresponding instructions for use of the implant manufacturer.

Note: DD Solid Connect® CAD/CAM Abutments of the respective series can only be combined with the corresponding compatible implant system. This also applies to the corresponding platform diameters.

Explanation of symbols:

	Manufacturer
	Date of manufacture
	Batch code
	Catalog number
	Keep dry
	Non-sterile
	Do not use if package is damaged
	Torque
	Consult electronic instructions for use www.dentaldirekt.de/IFU
	Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist
	Do not re-use
	Medical Device
	Unique Device Identification
	MR Conditional

Disposal:

Dispose of waste according to applicable legislation.

*Dental Direkt validated milling centers can be found at:
www.dentaldirekt.de/en/validated-milling-centers

Warranty:

Our products are subject to strict quality control and correspond to current technical standards. Please refer to our current warranty terms and conditions for further details.

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority to the Member State in which the user and/or patient is established.

Our products are subject to continual development, so we reserve the right to make changes. You will find the latest version of the instructions for use on our homepage at: www.dentaldirekt.de/en

This version replaces the previous.

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