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Instruction Manual

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Material-related characteristics

Chemical composition [wt. %]

Polymethylmethacrylate (PMMA)	100
-------------------------------	-----

Physical specifications

Flexural strength	[ISO 20795-2]	[MPa]	64
Flexural modulus	[ISO 20795-2]	[MPa]	1950
Fracture toughness	[ISO 20795-2]	[MPam ^{1/2}]	2,0
Charpy impact strength, +23°C	[ISO 179/1eU]	[kJ/m ²]	45
Density	[ISO 1183]	[g/cm ³]	1,17
Water sorption	[ISO 20795-2]	[µg/mm ³]	< 24
Solubility	[ISO 20795-2]	[µg/mm ³]	< 0,13
Residual monomer content	[ISO 20795-2]	[%]	< 0,7

Explanation of symbols

Manufacturer	Date of manufacture	Use-by date
LOT	REF	MD
Batch Code	Catalog number	Medical Device
Keep dry	Keep away from sunlight	
Content (Quantity)	Consult instructions for use	

R_xonly

Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist.

1. Indications

DD Bio Splint P HI are transparent dental milling blanks out of PMMA for the manufacture of splints, therapeutic splints and bite regulators for long-term application in the oral cavity up to 12 months.

2. Contraindications

- intolerance to the contained components
- crowns, bridges, model casting with clasps, bars, narrow sublingual bars and transversal connectors

3. Safety information

Please pay attention to the information in the material safety data sheet in its current version. Avoid the inhalation of dust particles during processing. Wear protection gloves, safety glasses and a dust mask to prevent irritations of eyes, skin and respiratory system.

4. Handling and storage

The medical device may only be processed by trained personnel. Only use clean, defect-free material from undamaged packaging and ensure that the blanks are stored at room temperature, in their original cartons, in a dry location and protected from light. Avoid contact with fluids.

5. Milling and further processing

Milling must be carried out with a milling strategy adapted to PMMA. Only use one-edged milling tools in the milling machine or wet grinding procedures!

Further processing:

Tungsten carbide burs like the one you use for soft reline materials are suitable for rough adaption. For the finalization common cross-toothed millers can be used. Polishing should be carried out with regard to a low contact pressure, otherwise the resulting heat might lead to smearing or melting.

6. Cleaning

- ultrasonic bath (max. 40 °C) with mild cleaning agent

Do not use the following methods for cleaning:

- steam blasting
- autoclave method (e.g. sterilization / disinfection)
- organic solvents and acids (e.g. acetone, hydrochloric acid)
- highly concentrated alcohols
- chemical denture cleaners
- abrasive cleaning agents

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7. Disposal

Remains of PMMA blanks can be disposed of in the regular household garbage.

We reserve the right to make changes because of the continuous development of our products. Please find the current version of the instruction manual on our website: www.dentaldirekt.de/en This version replaces all previous versions.