DDShadeConcept®

DE Gebrauchsanweisung, **EN** Instructions for Use, **FR** Notice d'utilisation, **IT** Istruzioni per l'uso,
ES Instrucciones de uso, **NL** Gebruiksaanwijzing, **EL** Οδηγίες χρήσης για το, **PT** Instruções de utilização,
SV Bruksanvisning, **ET** Kasutusjuhend, **DA** Brugsanvisning, **FI** Käyttöohje, **PL** Instrukcja używania,
CS Návod k použití skenovacího tělíška, **HU** Használati útmutató, **RO** Instrucțiuni de utilizare,
SK Návod na použitie, **BG** Инструкция за употреба, **HR** Upute za uporabu, **GA** Treoracha úsáide,
LV lietošanas pamācība, **LT** Naudojimo instrukcija, **SL** Navodila za uporabo, **MT** Struzzjonijiet għall-użu

DE	DEUTSCH	Gebrauchsanweisung	3
EN	ENGLISH	Instructions for Use	14
FR	FRANÇAIS	Notice d'utilisation	25
IT	ITALIANO	Istruzioni per l'uso	36
ES	ESPAÑOL	Instrucciones de uso	47
NL	NEDERLANDS	Gebruiksaanwijzing	58
EL	ΕΛΛΗΝΙΚΑ	Οδηγίες χρήσης για το	69
PT	PORTUGUÊS	Instruções de utilização	80
SV	SVENSK	Bruksanvisning	91
ET	EESTI	Kasutusjuhend	102
DA	DANSK	Brugsanvisning	113
FI	SUOMI	Käyttöohje	124
PL	POLSKI	Instrukcja używania	135
CS	ČEŠTINA	Návod k použití skenovacího tělíska	146
HU	MAGYA	Használati útmutató	157
RO	ROMÂNĂ	Instrucțiuni de utilizare	168
SK	SLOVENSKY	Návod na použitie	179
BG	БЪЛГАРСКИ	Инструкция за употреба	190
HR	HRVATSKI	Upute za uporabu	201
GA	GAEILGE	Treoracha úsáide	212
LV	LATVIEŠU VALODA	lietošanas pamācība	223
LT	LIETUVIŲ KALBA	Naudojimo instrukcija	234
SL	SLOVENŠČINA	Navodila za uporabo	245
MT	MALTI	Struzzjonijiet għall-użu	256

1. Produktbeschreibung

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1. Product description

The Dental Direkt water-based coloring liquids are suitable for complete or partial coloring of dental restorations made of Dental Direkt zirconium dioxide prior to final sintering.

Dental Direkt coloring liquids consist of the DD Basic Shade, DD Pro Shade C, DD Pro Shade Z und DD Art Elements versions.

To enable the dentist to provide each patient with customized and aesthetically pleasing restorations, Dental Direkt coloring liquids DD Basic Shade, DD Pro Shade C and DD Pro Shade Z are suitable for producing the basic tooth color and DD Art Elements for customizing the shade.

2. Intended purpose

Dental Direkt coloring liquids are intended for matching the color of dental restorations made of Dental Direkt zirconium dioxide milling blanks to the natural tooth color.

3. Indications

DD Basic Shade

For coloring zirconium dioxide restorations made of DD Bio Z, DD Bio ZX², DD cube ONE® or DD cubeX²® milling blanks.

DD Pro Shade Z

For coloring zirconium dioxide restorations made of DD Bio Z, DD Bio ZX² or DD cube ONE® milling blanks.

DD Pro Shade C

For coloring zirconium dioxide restorations made of DD cubeX²® milling blanks.

DD Art Elements

For coloring and customizing zirconium dioxide restorations made of DD Bio Z, DD Bio ZX², DD cube ONE® or DD cubeX²® milling blanks.

Zirconium dioxide restorations made of already precolored blanks must only be treated with DD Art Elements and not additionally with DD Basic Shade, DD Pro Shade C or DD Pro Shade Z.

4. Contraindications

Do not use in the event of intolerance to the constituents. Observe the contraindications of the Dental Direkt zirconium dioxide used in each case. The coloring liquids are not suitable for zirconium dioxide ceramics from other manufacturers.

5. Intended users

Dental Direkt coloring liquids must only be used by dentists and dental technicians in compliance with the procedures in the Instructions for Use.

6. Intended patient group

Coloring liquids for matching the shade of dental restorations made of Dental Direkt zirconium dioxide milling blanks are suitable for the permanent dentition of adult patients of any gender and nationality.

7. Handling and storage

Before using the Dental Direkt coloring liquid for the first time, check that the bottle is intact. Always store Dental Direkt coloring liquids in the tightly closed original container at a temperature no higher than 25°C, and ideally in the refrigerator (6-8°C). After using the DD Colorbox, the Dental Direkt coloring liquid must be poured back into the original container for storage.

To avoid contamination, Dental Direkt coloring liquid that has already been used should not be mixed with unused Dental Direkt coloring liquid. Failure to follow these rules may result in color changes.

8. Instructions for Use in the laboratory

Caution: Wear gloves, protective goggles and a face mask to avoid skin irritation. The liquid is not suitable for consumption. Please do not eat, drink or smoke when using the Dental Direkt coloring liquid.

The product should only be used by trained personnel and using proven techniques.

8.1 Application overview of DD Basic Shade, DD Pro Shade C and DD Pro Shade Z

For all three versions, 3 Bleach and 16 VITA® shades are available.

8.2 Application overview of DD Art Elements

For DD Art Elements, 10 effect and 6 incisal colors are available, plus the lightening liquid Vario.


Zirconium dioxide restorations made of already precolored blanks must only be treated with DD Art Elements and not additionally with DD Basic Shade, DD Pro Shade C or DD Pro Shade Z.

Incisal colors

There are 6 incisal colors for occlusal and incisal customization to choose from:

SA1	reddish / brownish shade:	for brush technique (OM3 · A1 · A2 · B1 · B2)
SA2	reddish / brownish shade:	for brush technique (A3 · A3.5 · A4 · B3 · B4)
SC1	grayish shade:	for brush technique (C1 · C2 · D2)
SC2	grayish shade:	for brush technique (C3 · C4 · D3 · D4)
S0	neutral shade:	for brush technique (OM1 · OM2) for immersion technique (all color versions)
S0 Bright	neutral shade / intense brightness:	for brush technique (OM1 · OM2) for immersion technique (all color versions)

The more intense brightness in the incisal area with the incisal color S0 Bright, compared to the S0 incisal color, is illustrated below:

* The zirconium dioxide specimens made of DD cube ONE® were immersed in the coloring liquid DD Pro Shade Z col. A4 for 15 seconds.

Effect colors

10 highly concentrated effect colors are available for tinting DD Basic Shade, DD Pro Shade C or DD Pro Shade Z.

The effect colors can also be used on their own to characterize fissures, incisal areas or gingival areas, for example.

When using DD Art Elements effect colors for mixing with DD Basic Shade, DD Pro Shade C or DD Pro Shade Z, a defined quantity should always be added to an entire packaging unit (30 ml or 100 ml) of DD Basic Shade, DD Pro Shade C or DD Pro Shade Z. This is the only way to ensure reliably reproducible shade results.

Vario: lightening liquid

The intensity of the coloring liquids is easily toned down to suit individual preferences by mixing with DD Art Elements Vario. For reproducible shade results, the dosing syringe included in the DD Art Elements Vario set is recommended for diluting 30 ml or 100 ml of DD Basic Shade, DD Pro Shade C or DD Pro Shade Z. The information below is intended as a basic guide for diluting DD Basic Shade, DD Pro Shade C or DD Pro Shade Z with DD Art Elements Vario. It is advisable to make color samples before using the mixtures.

Examples of mixing ratios for grading DD Basic Shade / DD Pro Shade C / DD Pro Shade Z in %						
DD Art Elements Vario	0 %	5 %	10 %	15 %	20 %	25 %
DD Basic Shade / DD Pro Shade C / DD Pro Shade Z	100 %	95 %	90 %	85 %	80 %	75 %

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8.3 Pigmentation of the coloring liquids

All the colors have a special pigmentation to make coloring easier for the user. This pigmentation has no effect on the shade that results after sintering.

DD Basic Shade, DD Pro Shade C and DD Pro Shade Z (all color versions):

Pigmentation in red.

DD Art Elements incisal colors (SA1, SA2, SC1, SC2, SO, SO Bright):

Pigmentation in yellow.

DD Art Elements effect colors: light pink, purple, yellow, orange

Pigmentation in orange.

DD Art Elements effect colors: light brown, dark brown

Pigmentation in green.

DD Art Elements effect colors: light grey, dark grey, graphite

Pigmentation in gray-violet.

DD Art Elements effect color: blue

Pigmentation in blue.

8.4 Coloring tools


Plastic tweezers

DD Bio Z applicator

Brush

DD Colorbox

Plastic instruments should be used to avoid metal marks.


Either the DD Bio Z applicator or a brush without a metal ferrule can be used for application.

To avoid contaminating the coloring liquid, keep the DD Colorbox closed when not in use.

The original container, which can be tightly closed, is recommended for storing the coloring liquid.

8.5 Immersion technique

If the restoration has been wet milled or wet ground (e.g. when using DD Bio ZS or DD Bio ZX² 19), it should be dried in a (ceramic) furnace to remove the cooling and/or lubricating liquid from the porous structure. Drying in the furnace is a prerequisite for a homogeneous shade result when using Dental Direkt coloring liquids.

DD Bio ZS – we recommend:

Drying the frames at 700°C (for 5 mins)


DD Bio ZX² 19 – we recommend:

Drying the frames at 80°C (for 30 mins) or 150°C (for 10 mins).

Shake the bottle well before use and wear nitrile gloves to avoid grease marks or contamination!

Immersion technique step-by-step

1-2x

1x 2x

DD Art Elements **SO Bright**

DD Art Elements **SO**

1. Applying Incisal SO / SO Bright


- The incisal color can be applied using a suitable applicator (e.g. DD Bio Z applicator) or a brush
- Make sure that the applicator is moderately wetted to prevent excess color liquid from running over the restoration
- In the upper cusp or incisal edge area, apply incisal color 1 to 2 times

Applying too much incisal color will lead to brightening of the DD Basic Shade, DD Pro Shade C or DD Pro Shade Z altogether. The dashed blue line on the crown shows approximately the furthest point to which the incisal color should be applied.

The restoration should be immersed as soon as the incisal color has been applied!


2. Preparing for immersion

- DD Colorbox is recommended for immersion

TIP:





When using less than 100 ml of coloring liquid, the necessary immersion volume can be achieved simply by tilting the DD Colorbox.

Recommended immersion time: 20 sec.

3. Immersion process

- The restoration is carefully placed in the sieve insert of the DD Colorbox using plastic tweezers
- Immersion time can be 15-60 sec. For a perfect result we recommend 20 sec.

 <p>  Incisal color  Incisal / dentin transition  DD Basic Shade, DD Pro Shade C or DD Pro Shade Z </p>	<ul style="list-style-type: none"> • When the restoration is immersed, the dentin and incisal colors are mixed, creating a transition between the two areas • Once the immersion time has elapsed, the restoration is simply lifted out of the DD Colorbox using the sieve insert • Alternatively plastic tweezers can be used
	<ul style="list-style-type: none"> • Carefully dab off excess coloring liquid (especially inside the crown and in the occlusal area) with an absorbent cloth, or blow it off with compressed air
<p> The colored restoration can be placed on a clean, absorbent surface (e.g. kitchen towel) for a short time to drip dry, but must not be stored on it.</p>	
<h2 style="background-color: orange; color: white; padding: 5px;">TIP 1</h2>	
	<p>Combining with the brush technique enables customization when using the immersion technique.</p> <ul style="list-style-type: none"> • Apply incisal color SO / SO Bright as described previously • After applying the incisal SO / SO Bright in the upper cusp or incisal edge area, the tooth neck can also be colored. • Depending on the existing tooth shade, the following options are available for coloring the tooth neck: DD Art Elements light brown / dark brown, DD Art Elements orange, DD Basic Shade / DD Pro Shade C / DD Pro Shade Z (one shade darker)
 <p>Pre-drying: 10 mins.</p> 	<ul style="list-style-type: none"> • Carry out immersion in DD Basic Shade, DD Pro Shade C or DD Pro Shade Z • Then pre-dry for 10 mins.
	<ul style="list-style-type: none"> • The pre-dried crown can now be further customized with effect colors in the “blocked incisal area” and is then pre-dried, as described under “Drying”, before sintering <p>TIP: Use the DD Art Elements effect colors for customization.</p>


TIP 2

To avoid extensive infiltration of DD Basic Shade, DD Pro Shade C or DD Pro Shade Z when using the immersion technique, the areas on solid pontics can be blocked with DD Art Elements incisal colors.

- With highly translucent zirconium dioxide in particular, blocking can prevent pontics from appearing darker than thin-walled crowns or copings after a long immersion time
- Before immersion, apply 1-2 thick drops of DD Art Elements SO / SO Bright basally and allow to infiltrate slowly

8.6 Brush technique

If the restoration has been wet milled or wet ground (e.g. when using DD Bio ZS or DD Bio ZX² 19), it should be dried in a (ceramic) furnace to remove the cooling and/or lubricating liquid from the porous structure. Drying in the furnace is a prerequisite for a homogeneous shade result when using Dental Direkt coloring liquids.

DD Bio ZS – we recommend:

Drying the frames at 700°C (for 5 mins)






DD Bio ZX² 19 – we recommend:

Drying the frames at 80°C (for 30 mins) or 150°C (for 10 mins.)


Shake the bottle well before use and wear nitrile gloves to avoid grease marks or contamination!

Brush technique step-by-step

 <p>1-2x 2-3x</p>	<p>1. Applying dentin color</p> <ul style="list-style-type: none"> • Application of DD Basic Shade, DD Pro Shade C or DD Pro Shade Z
 <p>1x</p>	<p>2. Applying incisal color (SA / SC colors)</p>
	<p>3. Customizing with effect colors</p> <ul style="list-style-type: none"> • 1-2 times dark brown / light brown (select according to tooth shade)
	<p>Sequence for applying effect colors:</p> <ul style="list-style-type: none">  1. Dark brown, light brown effect  2. Blue effect  3. Light grey, dark grey, graphite effect  4. Light pink, purple, yellow, orange effect <p>Accentuation with effect colors:</p> <p>Shaping (colored) fissures, crown margins</p> <ul style="list-style-type: none"> • Yellow, orange, light brown, dark brown effect <p>More depth for tooth cusps and incisal edges</p> <ul style="list-style-type: none"> • Purple, blue, light gray, dark gray, graphite effect <p>Shaping of gingival areas</p> <ul style="list-style-type: none"> • Light pink effect

Before using a new color, carefully wash out and dry the brush!


4. Applying dentin color to inner crown surfaces

To achieve the best shade result, the inner surface of the crown must finally be colored with DD Basic Shade, DD Pro Shade C or DD Pro Shade Z. The coloring liquid should be applied liberally several times in order to infiltrate the material completely.

Store the colored restoration on a non-absorbent surface!



Painting solid parts

For solid pontics, the frequency of color application should be increased. This also applies to crowns, partial crowns and caps with thick walls.

8.7 Drying

For pre-drying the restorations, a small household fan oven with a timer function is recommended, as this ensures virtually uniform heat distribution, in contrast to other drying methods.

Drying time: 10–30 Min.
Temperature: 80–100 °C

Place your frames on a ceramic surface in the furnace so that they can be “bathed” by the circulating air.

To achieve a homogeneous shade result, the frames should be dried straight after coloring. For this purpose, we recommend positioning crowns and bridges on the occlusal surface.

You should use this furnace exclusively for drying your zirconium dioxide restorations and not for preparing food!

8.8 Sintering

Please refer to our separate sintering instructions!

Sintering cycle with normal furnace filling, without cover:

- ↑ Heating to 900°C (8°C / min.),
- 30 mins hold time at 900°C,
- ↑ Heating up to end temp. 1450°C (3°C / min.),,
- 120 mins hold time at 1450 °C,
- ↓ Cooling to at least 200°C (10°C / min.)

8.9 Glazing

An additional glaze application not only results in an exact match with the desired shade, but also makes the restoration look more natural and aesthetic.

Firing charts for DD contrast® glaze and DD Nature ZR Glaze

Firing	Start temp. [°C]	Drying time [mins]	Closing time [mins]	Heating rate [°C/min.]	End temp. [°C]	Holding time [mins]	Cooling [mins]	Vacuum
Glaze firing	450	3	3	45	810	1	3	Max. (= 100%)
DD contrast® glaze								
Glaze firing	440	2	–	45	815	1	3	off
DD Nature ZR Glaze								
Full anatomy glaze firing	440	2	–	45	830	1	3	off

We recommend the use of glazes fired at a maximum of 850°C.

Increased temperatures as well as multiple repetitions of glaze firing can lead to color changes.

TIP 3

Color assessment on the model

- Light radiation affects color perception
- Copings are more translucent than the solid pontic, which can make the solid pontic appear darker

Always place the bridge on the model for color assessment!

9. Possible side effects and interactions

No known side effects or interactions.

10. Disposal

Dispose of residual quantities of the coloring liquid and the bottle in accordance with the applicable regulations. The coloring liquid must not be tipped down the drain or disposed of with household waste.

Please note: Please take note of the information in the latest version of the safety data sheet.

11. Reporting incidents

Any serious incident that occurs in connection with the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located.

Please note: The safety and clinical performance summary report can be requested at info@dentaldirekt.de.

We are continuously developing and enhancing our devices, and therefore reserve the right to make changes. The latest version of the Instructions for Use can be found on our website at:

www.dentaldirekt.de/en/IFU

This version replaces all previous versions.

Explanation of symbols:

Manufacturer

Date of manufacture

Use before

Batch code

Catalog number

Keep away from sunlight

Upper temperature limit

Consult electronic Instructions for Use
www.dentaldirekt.de/en/IFU

Caution: Under U.S. federal law, the device may only be sold by or on behalf of a dentist.

Medical device

Health hazard (see safety data sheet for more information)

Corrosive (see safety data sheet for more information)

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